
# **Rapid COVID19 IgG/IgM Test**

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
## BIBLIOGRAFIA/REFERENCES

1. Schmaljohn C, Hjelle B. Coronavirus 19 (COVID19) es: a global disease problem. *Emerging Infectious Diseases* 1997;3:95-104.
2. Khan A, Ksiazek TG, Peters CJ. Viral hemorrhagic fevers. *Seminars in Pediatric Infectious Diseases* 1997;8(Suppl 1):64-73.
3. Peters CJ. Viral hemorrhagic fevers. In: *Viral Pathogenesis*. New York City: Lippincott-Raven Publishers, 1997, 779-94.
4. Peters CJ, Simpson GL, Levy H. Spectrum of Coronavirus 19 (COVID19) infection: hemorrhagic fever with renal syndrome and Coronavirus 19 (COVID19) pulmonary syndrome. *Annual Review of Medicine* 1999;50:531-45.
5. Linderholm M, Elgh F. Clinical characteristics of Coronavirus 19 (COVID19) Infections on the Eurasian continent. *Current Topics in Microbiology and Immunology*. 2001;256:135-51.

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*Is compliant with the provisions of the Council Directive 98/79/CE for in vitro diagnostic medical devices. All supporting documentation is retained under the premises of the manufacturer.*

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## **SECTION 1: IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING**

### **1.1. Product identifier**

Applicable Product REF:

Product name: **Rapid New Coronavirus 19 (COVID19) IgG/IgM combo test**

Composition: Test + Sample buffer

### **1.2. Relevant identified uses of the substance or mixture and uses advised against**

*Relevant identified uses:*

- TEST DEVICE: Medical device for professional *in vitro* diagnostic use only. Use for detection of *Novel Coronavirus 2019 (COVID19)* in whole blood, serum and plasma samples, indicative of *Novel Coronavirus 2019 (COVID19)* infection.
- SAMPLE DILUENT: Use for Anti *Novel Coronavirus 2019 (COVID19)* IgG and IgM antibodies extraction contained in whole blood, serum and plasma samples. This buffer is only provided with the product which has to be used (**Rapid New Coronavirus 19 (COVID19) IgG/IgM combo test**).

*Uses advised against:* No information available.

### **1.3. Manufacturer identification**

Company:

Registered office:

Telephone: ; Fax:

E-mail:

For additional information:

### **1.4. Emergency line:**

## **SECTION 2: HAZARD IDENTIFICATION**

**2.1 Classification of the substance or mixture:** Non-hazardous preparation (Regulation 1272/2008/EC).

**2.1.1. Classification according to Regulation (EC) No 1272/2008 [CLP]:** Non-hazardous.

**2.1.2. Classification according to Directive 1999/45/EC:** Non-hazardous.

**2.1.3 Additional information:** See SECTION 16.

### **2.2 Label elements**

Signal Word: **None**

**2.3 Other hazards:** No hazards know.

## **SECTION 3: COMPOSITION / INFORMATION ON INGREDIENTS**

**3.1 Substances:** No information available.

### **3.2 Mixtures**

	<b>MSDS (Material Safety Data Sheet)</b>	<b>Rev. 0</b>
	Pag. 8 di 12	10/03/2020

**Mixture description:** Sample diluent contains buffer, salt and <0.1% of sodium azide as preservative.

### 3.2.1 Hazardous components:

Substance	CAS No.	CE No.	W/W %	Regulation (EC) no1272/2008 [CLP]	Classification
Sodium azide (Preservative)	26628-22-8	247-852-1	<0.1	Acute Tox. 2: H300 Aquatic Acute 1: H400 Aquatic Chronic 1: H410	Due to concentration <0.1%, this preparation is not classified as dangerous on the basis of health and/or environment effects.

Lowest generic cut-off value:  $\geq 0.1$

Lowest specific concentration limits/ M-factor: N/A (according to ATE Annex I section 3.1.3.6.1. and Table 3.1.2, classification  $\geq 1.0$  %).

**Note:** Sample buffer is classified as a not dangerous preparation (Regulation (EC) No 1272/2008 [CLP]).

The device consists in a strip composed of several layers: an absorbent material pre-dried with a coloured latex-antibodies conjugate against the product antibodies, a nitro-cellulose membrane with coated antigens against the product antibodies and cellulose absorbent. Contains sodium azide <0.1% as preservative.

**Additional information:** For full text of H-phrases: see SECTION 16.

## **SECTION 4: FIRST-AID MEASURES**

### **4.1 Description of first aid measures**

**Following eye contact:** Rinse thoroughly with plenty of water for at least 15 minutes. Consult a physician.

**Following skin contact:** Wash off immediately with soap and plenty of water. Consult a physician.

**Following ingestion:** Clean mouth with water and drink afterwards plenty of water. Consult a physician.

**Following inhalation:** Ensure sufficient ventilation of workplace. Consult a physician.

**4.2 Most important symptoms and effects, both acute and delayed:** No information available.

**4.3 Indication of any immediate medical attention and special treatment needed:** Treat symptomatically.

## **SECTION 5: FIREFIGHTING MEASURES**

### **5.1 Extinguishing media**

**Suitable Extinguishing Media:** Water or CO<sub>2</sub>. Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

**Extinguishing media which must not be used for safety reasons:** No information available.

**5.2 Special hazards arising from the substance or mixture:** Thermal decomposition can lead to release of irritating gases and vapors.

**5.3 Advice for firefighters:** As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

## **SECTION 6: ACCIDENTAL RELEASE MEASURES**

	<b>MSDS (Material Safety Data Sheet)</b>	<b>Rev. 0</b>
	<b>Pag. 9 di 12</b>	<b>10/03/2020</b>

**6.1 Personal precautions, protective equipment and emergency procedures:** Prevent contact with skin, eyes and clothes. Use personal protective equipment. Ensure adequate ventilation.

**6.2 Environmental precautions:** Given the way dispensation there is no possibility of accidental spillage in sufficient quantity to be dangerous. Avoid release to the environment.

**6.3 Methods and material for containment and cleaning up:** Soak up with inert absorbent material. Clean contaminated surface thoroughly.

**6.4 Reference to other sections:** If appropriate Sections 8 and 13 shall be referred to.

#### **SECTION 7: HANDLING AND STORAGE**

**7.1 Precautions for safe handling:** Good Laboratory Practices (disposal gloves). Not to eat, drink and smoke in work areas. Avoid contact and contamination with skin, eyes and clothes. Use disposal gloves. Specimens should be handled as potentially infectious materials.

**7.2 Conditions for safe storage, including any incompatibilities:** Store in a dry place at +2°C to +30°C. Avoid storage near to heat sources.

**7.3 Specific end use(s):** Only use provided buffer for sample dilution.

#### **SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION**

**8.1 Control parameters:** Any specific protection and prevention measures should not be taken during use of the product.

##### **Exposure limits:**

<b>Substance</b>	<b>LTEL (8 hr)</b>	<b>STEL</b>
Sodium azide (NaN <sub>3</sub> )	<b>MG/m<sup>2</sup></b>	<b>MG/m<sup>2</sup></b>
<b>CAS No. 26628-22-8</b>	0.1	0.3

**8.2 Exposure controls:** All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.

**8.2.1 Appropriate engineering controls:** No relevant for this material.

**8.2.2 Personal protective equipment:** Handle with disposable gloves (EN 374). Wear appropriate protective safety eyewear and clothing, such as a lab coat.

**8.2.3 Environmental Exposure Controls:** No special measures are required.

#### **SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES**

##### **9.1 Information on basic physical and chemical properties**

###### **Appearance/Physical State**

**Test:** White solid with yellow lines. Solid reaction strip that is packaged depending on the format used.

**Sample diluent:** Transparent slightly yellowish.

	<b>MSDS (Material Safety Data Sheet)</b>	<b>Rev. 0</b>
	Pag. 10 di 12	10/03/2020

The following table only applies to sample diluent:

<b>Odor</b>	Odourless	<b>Explosion Limits</b>	Not applicable
<b>pH</b>	8.0-9.0	<b>Vapor Density</b>	Not determined
<b>Boiling Point</b>	Similar to water (100°C)	<b>Relative density</b>	Similar to water (1g/cm <sup>3</sup> )
<b>Flash Point</b>	Not applicable	<b>Solubility</b>	Soluble
<b>Vapor Pressure</b>	Similar to water (23hPa)	<b>Flammability</b>	Not applicable
<b>Melting Point</b>	Similar to water (0°C)	<b>Viscosity</b>	Not determined
<b>Autoignition Temp.</b>	Not determined	<b>Explosive Properties</b>	Not explosive
<b>Partition Coefficient (n-octanol/water)</b>	Not determined	<b>Oxidizing Properties</b>	Not determined

### **SECTION 10: STABILITY AND REACTIVITY**

**10.1 Reactivity:** No hazardous reactivity known.

**10.2 Chemical stability:** Under storage at normal ambient temperatures the product is stable. No known hazardous reactions.

**10.3 Possibility of hazardous reactions:** Thermal decomposition can lead to release of irritating gases and vapors.

**10.4 Conditions to avoid:** Direct contact with a flame. Temperatures outside the range of 2-30 ° C. Avoid storing in places with high humidity.

**10.5 Incompatible materials:** The sample should be treated only with buffer that is provided with the product before testing.

**10.6 Hazardous decomposition products:** No known hazardous decomposition products.

### **SECTION 11: TOXICOLOGICAL INFORMATION**

#### **11.1 Information on toxicological effects**

**Acute toxicity:** Product does not present an acute toxicity hazard based on known or supplied information. Oral LD50 Rat: 27mg/kg; Dermal LD50 Rabbit: 20mg/kg.

**Skin corrosion/irritation:** Based upon the available data, the classification criteria are not met.

**Serious eye damage/irritation:** Based upon the available data, the classification criteria are not met.

**Respiratory or skin sensitisation:** Based upon the available data, the classification criteria are not met.

**Germ cell mutagenicity:** Based upon the available data, the classification criteria are not met.

**Carcinogenicity:** A4-Not classifiable as a Human Carcinogen.

**Reproductive toxicity:** Based upon the available data, the classification criteria are not met.

**Summary of evaluation of the CMR properties:** Based upon the available data, the classification criteria are not met.

**STOT-single exposure:** Based upon the available data, the classification criteria are not met.

**STOT-repeated exposure:** Based upon the available data, the classification criteria are not met.

**Aspiration hazard:** Based upon the available data, the classification criteria are not met.

	<b>MSDS (Material Safety Data Sheet)</b>	<b>Rev. 0</b>
	<b>Pag. 11 di 12</b>	<b>10/03/2020</b>

## **SECTION 12: ECOLOGICAL INFORMATION**

**12.1 Toxicity:** Based upon the available data, the classification criteria are not met. The product should be discarded in a proper biohazard container after testing. Do not allow product to reach ground water, water bodies or sewage system.

**12.2 Persistence and degradability:** Based upon the available data, the classification criteria are not met.

**12.3 Bioaccumulative potential:** Based upon the available data, the classification criteria are not met.

**12.4 Mobility in soil:** Based upon the available data, the classification criteria are not met.

**12.5 Results of PBT and vPvB assessment:** No data available for assessment.

**12.6 Other adverse effects:** Based upon the available data, the classification criteria are not met.

## **SECTION 13: DISPOSAL CONSIDERATIONS**

### **13.1 Waste treatment methods**

**Waste from Residues:** After testing, the product must be disposed of compliance with the respective local, state or national regulations. One option would be possible inactivation of infectious agents in the product after use. Performed in autoclave at a pressure and a certain temperature.

**Non-contaminated packaging:** The containers can be recycled.

## **SECTION 14: TRANSPORT INFORMATION**

**Maritime transport (IMDG/IMO):** Not dangerous preparations not required transport regulations.

**Land transport (ADR):** Not dangerous preparations not required transport regulations.

**Air Transport (IATA):** Not dangerous preparations not required transport regulations.

## **SECTION 15: REGULATORY INFORMATION**

**15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture:** This product does not require special labelling, in accordance with the appropriate EC directives. These products are used for *in vitro* diagnosis, so they must meet the criteria described in Directive 98/79/CE, do not carry the CE marking for marketing outside the EU.

The product is a mixture which is not subject to Regulation (EC) No 1005/2009, (EC) No 850/2004.

**National Regulations:** Please ask your national/regional authorities.

**15.2 Chemical Safety Assessment:** A Chemical Safety Assessment/Report has not been conducted.

## **SECTION 16: OTHER INFORMATION**

**Recommendations:** Consult instructions for use prior to product use. Professional use only for *in vitro* diagnosis.

	<b>MSDS (Material Safety Data Sheet)</b>	<b>Rev. 0</b>
	<b>Pag. 12 di 12</b>	<b>10/03/2020</b>

**References (previous version):** approving the Regulation on classification, packaging and labeling of dangerous preparations, which incorporates into Italian law Directive 1999/45/CE, Directive 2001/60/CE and partly Directive 2001/58/CE. Directive 91/155/CE.

**Changes:** Update in accordance with Regulation (EC) No 1272/2008 and EU No 2015/830 (changes to all sections).

**Abbreviations and acronyms:**

**STOT:** Specific Target Organ Toxicity **Acute Tox. 2:** Acute toxicity: Category 2

**STEL:** Short Term Exposure Limit **Aquatic Acute 1:** Hazardous to the aquatic environment Acute: Category 1

**LTEL:** Long Term Exposure limit **Aquatic Chronic 1:** Chronic Aquatic toxicity: Category 1

**GHS:** Global Hazard Symbol

**PBT:** Persistent, Bioaccumulative and Toxic

**vPvB:** very Persistent and very Bioaccumulative

**ATE:** Acute toxicity estimate

**LD** Lethal Dose

**Key literature references and sources for data:** see instruction for use, Safety Data sheet and ECHA.

**Classification and procedure used to derive the classification for mixtures according to Regulation (EC) 1272/2008 [CLP]:** Annex I section 3 and 4; Annex VI Table 3.1 of Regulation (EC) No 1272/2008 was used for the purpose of classification.

**Relevant H-statements (number and full text):**

**H300:** Fatal if swallowed

**H400:** Very toxic to aquatic life.

**H410:** Very toxic to aquatic life with long lasting effects

**Training advice:** No special training is required.

**Contacts**

Via

Tel.

E-mail:

For more information visit:

**HEREIN BELIEVES THE INFORMATION PRESENTED IS ACTUAL AS IT HAS BEEN DRIVEN FROM THE WORKS AND OPINIONS OF PERSONS TO BE BELIEVED TO BE QUALIFIED EXPERTS.**

**HOWEVER, NOTHING CONTAINED IN THIS INFORMATION IS TO BE TAKEN AS A WARRANTY OR REPRESENTATION WHICH BEARS LEGAL RESPONSIBILITY. THE USER SHOULD REVIEW ANY RECOMMENDATIONS IN THE SPECIFIC CONTEXT OF THE INTENDED USE TO DETERMINE WHETHER THEY ARE APPROPRIATE.**

# COVID-19

C  
IgG  
IgM

C  
IgG  
IgM

S

One step rapid test

# COVID-19 IgG/IgM

## Rapid Test Cassette

(Whole Blood/Serum/Plasma)

Contents: Test Cassettes - 30  
Droppers - 30

Package Insert - 1  
Buffer - 2

[REF]

[LOT]

2020-03

2022-03

IVD

CE

QC Passed